


Stress & Recovery Analyzer Runs On Mac

BRIEF INSTRUCTION TO INSTALL THE SRA SOFTWARE ON A MAC

Introduction

- ❑ The Stress And Recovery Analyzer (SRA) is a native Windows application.
- ❑ SRA operates on PCs running Windows 10 or 11.
- ❑ It is possible to run the BHA on Mac computers with a little investment of money and time to set it up.
- ❑ The solution is to use a Mac program that creates a virtual PC inside a Mac.

Virtual PC




- ❑ Your Mac PC must run a special app that creates a Virtual PC.
- ❑ Virtual PC runs Windows 10 or 11 (Guest OS).
- ❑ Windows runs the BHA software.

Parallels Desktop for Mac



- ❑ Parallels Desktop is the native Mac software
- ❑ It creates a Virtual PC imitating a real PC device inside Mac.
- ❑ It runs Windows 10/11 inside the virtual PC.
- ❑ Any Windows software may run inside Windows on Mac.

Parallels Is Not Free



- ❑ Parallels is a subscription service
- ❑ Parallels costs \$100 per year.
- ❑ Windows 10/11 license must be purchased to run in Parallels.
- ❑ 14-days trial is available for Parallels and Windows license.

Download, Install, and Set Up Parallels




[Watch this video](#) to learn how to set up Parallels and Windows 11 on your Mac.

Parallels Running on Mac




- ❑ When Parallels Desktop with Windows is installed, it could be already running.
- ❑ Check the dock bar at the bottom to see if Parallels and Windows already running.
- ❑ Click on it to bring it forward.

Parallels Not Running on Mac




- ❑ If Parallels is not currently running, double-click on its shortcut icon located on the desktop screen.
- ❑ It will bring the Mac window containing Windows OS forward and start Windows.

Download SRA




- Go to your mailbox and open an order confirmation email.
- Click on the web link in it to start software downloading.

Download SRA




- ❑ Remember the location where the software installer will be downloaded.
- ❑ Click Save.

Download SRA




- ❑ Check if this icon appears in the Mac Downloads folder.

Install SRA




- ❑ Go to the Mac window containing running Windows.
- ❑ Click on the Desktop icon located on the task bar.

Install SRA




- Wait for the Home window to open in Windows.
- Click Downloads to open the Downloads folder.

Install SRA




- Find the **BIASCLOUDInstall** file.
- Double-click on it to launch the Stress & Recovery Analyzer installer.

Install SRA




- ❑ Check for this window to open.
- ❑ It may open behind the topmost window. Close or minimize it to reveal this one.
- ❑ Click **Yes**.

Install SRA




Click **OK** to confirm language selection.

Install SRA




Click **Install** to start the Stress & Recovery Analyzer installation.

Install SRA




- ❑ Click **Finish** to complete installation.
- ❑ Once installation is done, the Stress & Recovery Analyzer will launch.

Register SRA Account




- Now you need to register with the Stress & Recovery Analyzer as a user.
- Click **Register**.


Register SRA Account



- Now you need to register with the Stress & Recovery Analyzer as a user.
- Click **Register**.

Register SRA Account



 Click I am a home user.

Register SRA Account

Windows 11

BODY HEALTH ANALYZER

Create user account

Please fill the form below then click CREATE.

User name* Contact phone

First name* Contact e-mail*

Last name* Contact address

Gender* Male Female

Date of Birth* Month Day Year Health notes


Enter password*

Re-enter password*

* - Mandatory data entries


- Fill out the registration form.
- Make sure all mandatory data field are filled out.
- Click **Create**.

Register SRA Account




 Click **OK**.

Register SRA Account




- Open your mailbox registered in your newly created account.
- Open a new confirmation email.
- Click on the web link in it to confirm registration.

Register SRA Account




- Get your new account registration confirmation message.
- Now go back to the Stress & Recovery Analyzer window.

Log into SRA Account



- ❑ Enter your new account username and password.
- ❑ Click **Login**.

Log into SRA Account



📄 Congratulations! Now you are ready to use the Stress And Recovery Analyzer.